jota

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
SILKY SMOOTH  
RESULTS EVEN IN THE HARDEST CASES

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JOTA AG is the manufacturer and full-service supplier of a wide range of rotary instruments. For over a century, we have been developing and producing high-quality instruments that are used successfully in more than 80 countries worldwide.

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We know how to adapt and make a point of understanding different cultures, applications and special features; as far as possible, our aim is to match our services to the customer's specific needs.

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We are pleased to present our podiatry catalogue. We have a specially selected range of rotary instruments which are perfect for use in the foot care sector.

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



The total lengths of instruments can be longer or shorter according to type of construction.

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## Maximum speed:

The maximum permissible speed is shown on every JOTA blister pack and must not be exceeded of safety reasons.

## Reference speed:



The speed recommendations given here are general reference speeds for use in the micromotor. Specifically recommended speeds can be found in the product-related information. The speed recommendations referred to must be adjusted to the particular applications and medical circumstances. Lower speeds should be selected when using the instruments dry to reduce any potential build-up of heat.

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These abrasive instruments are perfect for removing and smoothing calloused skin and nails. They feature a long service life, very good abrasive performance and a variety of possible uses. Careful manufacturing processes guarantee that the working parts, including the tips, radii and edges, are fully and densely coated. The natural diamonds used by JOTA are embedded using a galvanic process specially designed for use in podiatry. Their naturally very irregular geometries make them particularly abrasive. Diamond is the hardest material known and ensures that the instruments have a long service life. A high-quality, corrosion-resistant steel is used for the shanks, making the instruments suitable for all familiar sterilisation and cleaning methods. For intermediate cleaning, we recommend our diamond cleaning stone with which any blocked spaces between the individual diamond particles can be opened up again. CAUTION: This cleaning does not replace final reprocessing and sterilisation.

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The use of coarse-grain diamonds (ISO 533, 543 and 553) may lead to an increased build-up of heat. We particularly recommend irrigating with cooling liquid when using these instruments. When using the instruments dry, we recommend reducing the speed due to the build-up of heat and also to minimise the risk when treating high-risk patients.

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The perfect addition for skin removal, grinding and smoothing:

- › round edge on pass from top to side
- › safe treatment
- › perfect for scarf skin and nail surfaces

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Ideal for a safe therapy at nail borders and close to nail beds:

- › smooth, no-grinding top
- › special rounded-edge pass from top to side
- › safe treatment at border regions to skin and nail bed
- › no undesired skin contact
- › reduced risk for injuries
- › ideal for diabetics patients

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For a selective treatment of diseased, thickened or deformed nails:

- > on top diamond coating only
- > special rounded-edge pass from top to side
- > no contact of grinding surface with adjacent skin
- > reduced risk for injuries
- > ideal for diabetics patients

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Safe-Edge: round pass without sharp edges

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DiaKap: The permanent alternative for top-filing

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› to clear and remove persistent callosity and nail material  
› diamond coated with stronger, sharper edged grit  
› new structure preventing dogging effect of the sharp edged grit

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The impressive features of JOTA carbide instruments include their precise efficient cutting blades in various degrees of fineness, optimal concentric accuracy and long service life. The shank of the carbide cutters is made of stainless steel to ensure that they can be reprocessed easily. These cutting instruments are used for extensive removal and shaping of thickened or mycotic nails.

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JOTA steel instruments are made of stainless steel so they can be readily sterilised. The instruments posses a high degree of elasticity together with optimal material strength. The range comprises the most varied forms and blade designs, optimally suitable for removing thick layers of calloused skin, treating nail plates and corn removal.

#### PLEASE NOTE:

When sterilising stainless steel instruments, it should be noted that the instruments lose strength at temperatures above 180°C, significantly reducing their service life. JOTA AG recommends sterilising steel instruments only in an autoclave !

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Rose drill for work on corns, nails and nail-bed

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Rose drill with cross cut for work on corns, nails and nail-bed

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Fissure cutter with cross cut for working on nails and nail bed

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

specially for Mycose-Nails  
fast and gentle operation  
for thin out the nail-plate

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



### Special cutter:

- › hexagonal working part without sharp edges
- › specially suitable for gentle removal of deposits in the nail fold, cuticle and hangnails
- › does not leave any furrows behind, unlike a diamond instrument

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These abrasives, consisting of ceramic-bonded corundum, are mounted on a rustproof shank. This allows them to be cleaned and sterilised without difficulty. JOTA abrasives are optimally suitable for removing and contouring nail plates, nail edges and skin. The abrasives are designed to be used dry. Depending on the degree of soiling, the abrasives should be cleaned with a brush or in an ultrasonic bath. For heavy soiling, the **529 cleaning stone** can also be used beforehand.

Since effective sterilisation of the porous ceramic abrasives cannot always be guaranteed, we recommend using these instruments once only. For this reason, ceramic abrasives are not marketed as medical devices and therefore do not carry the CE mark !

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Abrasive caps are suitable for rapid removal of calloused skin by means of a suitable dry suction technique. Since the abrasive caps cannot be reprocessed or sterilised, they should be used exclusively as disposable articles. When working with abrasive caps, it is important to use the correct speed as otherwise excessive heat may be produced.

**Abrasive cups are designed for single use** and are not marketed by JOTA AG as medical devices. They therefore do not carry the CE mark !

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
The background of the entire image is a dense, irregular pile of small, yellowish-green stones or pebbles, possibly made of a mineral like chalcedony or jasper. They are piled high, filling the frame.

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JOTA polishers consist of silicon carbide abrasive grit bonded in an elastic silicon matrix. The shanks are made of stainless steel, allowing the instruments to be reprocessed.

The polishers are used to smooth and finish natural and artificial nails and nail prostheses and for specific finishing of nail edges. A smooth matte surface is obtained with the green and grey polishers. The beige produces a shiny surface and is therefore optimally suitable as a final step in manicures and pedicures. The goat hair brush is used for optimal cleaning of the nails and cuticles following treatment. A high-sheen surface can be achieved with the microfibre buffering wheel in combination with a polishing paste. The polishers, brushes and buffering wheels are designed for use in the dry technique.

We recommend a speed of 7'000 to 10'000 rpm.

Since polishers, brushes and buffering wheels are regarded by JOTA AG as „cosmetic products“, they do not carry the CE-mark.

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
Goat hair soft for polishing with polishing paste

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
Cleaning stone for diamonds  
(pre-cleaning of high soiled instruments)  
Dimensions: 75 x 25 x 8 mm

Moisten the stone slightly and draw the diamond instrument over it at a low speed. This does not replace sterilisation.

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# MANUAL AND MACHINE CLEANING PROCEDURE

Newly delivered unsterile instruments must be prepared accordingly before first use !

<b>1. Pre-treatment</b>	<ul style="list-style-type: none"> <li>Abrasive impurities need to be removed from the products directly after use (within two hours maximum). To do so, use running water or a disinfectant solution (must not contain aldehydes, its effectiveness should be established, suitable and compatible for the products)</li> </ul>	
<b>2a. Cleaning and disinfecting</b>	<ul style="list-style-type: none"> <li>Disassemble the products</li> <li>Place the disassembled products in the cleaning bath for the prescribed contact time and concentration, take care that the products do not touch each other</li> <li>Products with lumen: rinse all instrument lumens 5 times at the beginning and or at the end of the contact time using a disposal syringe and a cannula</li> <li>Take the products out of the cleaning bath and rinse at least 3 times thoroughly with water</li> <li>Check the products</li> <li>Dry the products by blowing them dry using filtered pressurised air (Use only filtered air for drying and water that is either sterile or low in germs and endotoxins)</li> <li>Wrap the products</li> </ul>	
<b>2b. Cleaning and disinfecting with disinfector / RDG (MACHINE CLEANING)</b>	<ul style="list-style-type: none"> <li>The effectiveness of the disinfecter has to be certified</li> <li>Disassemble the products as far as possible</li> <li>Place the disassembled products in the disinfecter and take care that the products do not touch each other</li> <li>Start the program</li> <li>Remove the products from the disinfecter after the program has finished</li> <li>Check and wrap the products straight after removal if possible</li> <li>Ensure that the disinfecter is regularly maintained and checked. Use only filtered air for drying and water that is either sterile or low in germs and endotoxins</li> </ul>	
<b>3. Checking / Maintenance / Packaging</b>	<ul style="list-style-type: none"> <li>Check the products for corrosion, damaged surfaces, bare patches, broken, chipped-off edges, deformations or impurities and eliminate damaged products. Products that are still contaminated need to be cleaned and disinfected once more</li> <li>Re-assemble disassembled products, instrument oil must not be used</li> <li>Arrange the cleaned and disinfected products in the dedicated bur block / sterilisation tray.</li> <li>Wrap the products, bur blocks, sterilisation trays using disposable sterilisation packaging (disposable or double packaging) or sterilisation containers</li> </ul>	
<b>4. Sterilisation</b>	<ul style="list-style-type: none"> <li>Steam sterilisation: use of a fractional vacuum process or a gravitation process</li> <li>Maximum sterilisation temperature 134°C (273°F)</li> <li>At least 3 min (or 18 min at prion deactivation) at 132°C (270°F) / 134°C (273°F) and drying time of 20 min</li> </ul>	
<b>5. Storage</b>	<ul style="list-style-type: none"> <li>The products need to be stored in sterilisation wrapping in a dry and dust-free place. Please note the storage-life resulting from the validation of the sterilisation wrapping</li> </ul>	
<b>i</b>	<ul style="list-style-type: none"> <li>New unsterile supplied instruments must be prepared before first use.</li> <li>Instruments made of tool steel are not suitable for sterilisation without an appropriate pretreatment.</li> <li>When cleaning instruments unprotected against corrosion (like steel instruments) a corrosion-inhibiting disinfectants and cleaning agents must be used and they have to be pretreated with rust preventing spray before sterilisation.</li> <li>Furthermore, make sure that instruments of different materials should never be reprocessed together.</li> <li>Use particular care during the cleaning of the grinding surfaces of the instruments and ensure that all residues are removed by using a brush.</li> <li>In case of severe contamination of the instrument, it is recommended to use an ultrasonic bath for cleaning.</li> <li>Usage of protection gloves during work with contaminated instruments is highly recommended!</li> <li>Bur blocks / instrument trays: <ul style="list-style-type: none"> <li>Cleaning and disinfecting only without products being loaded (products must not be cleaned and disinfected whilst they are in the bur block / instrument tray)</li> <li>Burblocks and other instruments made of aluminum (or other light metal) are not suitable for cleaning using a disinfecter RDG</li> </ul> </li> <li>Please follow the manufacturer's instructions and the maintenance specifications for using of the autoclave and the disinfecter.</li> <li>The method of use, reaction time and suitability of disinfectants and cleaning agents for certain types of instrument are covered by the manufacturers' instructions.</li> </ul>	
<p>These are a comprehensive description of our detailed reprocessing advices. This can be found also in other languages on our website, in our catalogue and on request we would be pleased to send you the printed information.</p>		

## Instructions for processing & General application and safety instructions

### **Instructions for the processing (cleaning, disinfection, and sterilization) of instruments from JOTA AG**

Issued: January 2017

The medical devices produced and sold by JOTA AG are re-usable unless their label contains explicit information to the contrary. However, as a rule, it is the sole responsibility of the doctor/expert using the devices to decide whether, depending on the respective case and the potential wear and tear of the products, he can re-use the products and how frequently he uses them. In case of doubt, it is always advisable to discard the products early and to replace them. The manufacturer JOTA AG cannot guarantee the faultless function and performance of the products combined with a maximum degree of safety if the products are overused. These reprocessing instructions apply in principle to all medical devices making up the product range supplied by JOTA AG. Any particular features and/or exclusions that only concern individual items or groups of items are referred to separately.

#### **Fundamental points**

All instruments are to be cleaned, disinfected, and sterilized prior to each application; this is required as well for the first use after delivery of the unsterile instruments (cleaning and dis-infection after removal of the protective packaging, sterilization after packaging). An effective cleaning and disinfection is an indispensable requirement for an effective sterilization of the instruments.

You are responsible for the sterility of the instruments. Therefore, please ensure that only sufficiently device and product specifically validated procedures will be used for cleaning, disinfection, and sterilization, that the used devices (WD, sterilizer) will be maintained and checked regularly, as well as that the validated parameters will be applied for each cycle.

Please pay attention to avoid a higher contamination of the complete bur block during application; otherwise it is necessary to clean and disinfect the bur block as well as all instruments inside (after removal).

Additionally, please pay attention to the legal provisions valid for your country as well as to the hygienic instructions of the doctor's practice or of the hospital. This applies particularly to the different guidelines regarding the inactivation of prions (not relevant for USA).

**Some instruments require additional aspects. For this, pay attention to chapter "Specific aspects".**

#### **Cleaning and disinfecting**

##### **Basic:**

If possible, an automated procedure (WD (Washer-Disinfector)) should be used for cleaning and disinfection of the instruments. A manual procedure – even in case of application of an ultrasonic bath – should only be used if an automated procedure is not available; in this case, the significantly lower efficiency and reproducibility of a manual procedure has to be considered.

The pre-treatment step is to be performed in both cases.

##### **Pre-treatment:**

Please remove coarse impurities of the instruments directly after application (within a maximum of 2 h).

##### **Procedure:**

1. Rinse the instruments at least 1 min under running water (temperature < 35 °C/95 °F).
2. Soak the instruments at least for the given soaking time in the pre-cleaning solution<sup>1</sup> (by the use of an ultrasonic bath) so that the instruments are sufficiently covered. Pay attention that there is no contact between the instruments. Assist cleaning by careful brushing with a soft brush (at least three times after beginning of soaking, aids see chapter "Specific aspects").
3. Activate ultrasonic treatment for an additional soaking time (but not less than 5 min).
4. Then, remove the instruments of the pre-cleaning solution and post-rinse them at least three times intensively (at least 1 min) with water.
5. In case of still visible contamination repeat steps 2, 3, and 4, otherwise discard the instrument. This is especially relevant for diamond instruments.

Pay attention to following points during selection of the cleaning detergent<sup>1</sup>:

- › fundamental suitability for the cleaning of instruments made of metallic or plastic material
- › suitability of the cleaning detergent for ultrasonic cleaning (no foam development)
- › compatibility of the cleaning detergent with the instruments (see chapter „material resistance“)

Pay attention to the instructions of the detergent manufacturer regarding concentration, temperature and soaking time as well as post-rinsing. Please use only freshly prepared solutions as well as only demineralized sterile or low contaminated water (max. 10 germs/ml) as well as low endotoxin contaminated water (max. 0.25 endotoxin units/ml), for example purified/highly purified water, and a soft, clean, and lint-free cloth and/or filtered air for drying, respectively.

<sup>1</sup> In case of application of a cleaning and disinfection detergent for this (e.g. in consequence of personnel's safety) please consider, that this should be aldehyde-free (otherwise fixation of blood impurities), possess a fundamentally approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking), be suitable for the disinfection of instruments made of metallic or plastic material, and be compatible with the instruments (see chapter „material resistance“). Please consider, that a disinfectant used in the pre-treatment step serves only the personnel's safety, but cannot replace the disinfection step later to be performed after cleaning.

#### **Automated cleaning/disinfection (WD (Washer-Disinfector)):**

Pay attention to following points during selection of the WD:

- › fundamentally approved efficiency of the WD (for example CE marking according to EN ISO 15883 or DGHM or FDA approval/clearance/registration)
- › possibility for an approved program for thermal disinfection (AO value > 3000 or – in case of older devices - at least 5 min at 90°C/194 °F; in case of chemical disinfection danger of remnants of the disinfectant on the instruments)
- › fundamental suitability of the program for instruments as well as sufficient rinsing steps in the program
- › post-rinsing only with demineralized sterile or low contaminated water (max. 10 germs/ml, max. 0.25 endotoxin units/ml), for example purified/highly purified water
- › only use of filtered air (oil-free, low contamination with microorganisms and particles) for drying
- › regularly maintenance and check/calibration of the WD

Pay attention to following points during selection of the cleaning detergent:

- › fundamental suitability for the cleaning of instruments made of metallic or plastic material
- › additional application – in case of non-application of a thermal disinfection – of a suitable disinfectant with approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) compatible to the used cleaning detergent
- › compatibility of the used detergents with the instruments (see chapter „material resistance“.)

Pay attention to the instructions of the detergent manufacturers regarding concentration, temperature and soaking time as well as post-rinsing.

##### **Procedure:**

1. Transfer the instruments in the WD by the use of a small pieces basket.
2. Start the program.
3. Remove the instruments of the WD after end of the program.
4. Check and pack the instruments immediately after the removal (see chapters „check“, „maintenance“, and „packaging“, if necessary after additional post-drying at a clean place).

The fundamental suitability of the instruments for an effective automated cleaning and disinfection was demonstrated by an independent, governmentally accredited and recognized (§ 15 (5) MPG) test laboratory by application of the WD G 7836 CD, Miele & Cie. GmbH & Co., Gütersloh, (thermal disinfection) and the pre-cleaning and cleaning detergent Neodisher mediclean forte (Dr. Weigert GmbH & Co. KG, Hamburg) considering to the specified procedure.

#### **Manual cleaning and disinfection:**

Pay attention to following points during selection of the cleaning and disinfection detergents:

- › fundamental suitability for the cleaning and disinfection of instruments made of metallic or plastic material
- › in case of application of an ultrasonic bath: suitability of the cleaning detergent for ultrasonic cleaning (no foam development)
- › application of a disinfectant with approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) compatible with the used cleaning detergent
- › compatibility of the used detergents with the instruments (see chapter „material resistance“.)

Combined cleaning/disinfection detergents should not be used.

Only in case of extremely low contamination (no visible impurities) combined cleaning/disinfection could be used.

Pay attention to the instructions of the detergent manufacturers regarding concentration, temperature and soaking time as well as post-rinsing. Please use only freshly prepared solutions as well as only demineralized sterile or low contaminated water (max. 10 germs/ml) as well as low endotoxin contaminated water (max. 0.25 endotoxin units/ml), for example purified/highly purified water, and a soft, clean, and lint-free cloth and/or filtered air for drying, respectively.

##### **Procedure:**

##### **Cleaning**

1. Soak the instruments for the given soaking time in the cleaning solution (by the use of a ultrasonic bath) so that the instruments are sufficiently covered. Pay attention that there is no contact between the instruments. Assist cleaning by careful brushing with a soft brush (at least three times after beginning of soaking, aids see chapter "Specific aspects").
2. Activate ultrasonic treatment for an additional soaking time (but not less than 15 min).
3. Then, remove the instruments of the cleaning solution and post-rinse them at least three times intensively (at least 1 min) with water.
4. Check the instruments (see chapters „check“, and „maintenance“).

##### **Disinfection**

5. Soak the instruments for the given soaking time in the disinfectant solution so that the instruments are sufficiently covered. Pay attention that there is no contact between the instruments.
6. Then, remove the instruments of the disinfectant solution and post-rinse them at least five times intensively (at least 1 min) with water.
7. Dry and pack the instruments immediately after the removal (see chapter „packaging“, if necessary after additional post-drying at a clean place).

The fundamental suitability of the instruments for an effective cleaning and disinfection was demonstrated by an independent, governmentally accredited and recognized (§ 15 (5) MPG) test laboratory by application of the pre-cleaning and cleaning detergent Cidzyme/Enzol and the disinfectant Cidex OPA (Johnson & Johnson

GmbH, Norderstedt) considering the specified procedure.

##### **Check**

Check all instruments after cleaning or cleaning/disinfection, respectively, on corrosion, damaged surfaces, and impurities. Do not further use damaged instruments (for limitation of the numbers of re-use cycles see chapter „reusability“). Still dirty instruments are to be cleaned and disinfected again.

##### **Maintenance**

Instrument oils or grease must not be use with the exception of steel instruments. In that case use only instrument oils (white oil) admitted to steam sterilization considering the maximum possible sterilization temperature, with approved biocompatibility and without mono-, di-, or triethanolamine as corrosion inhibitor.

##### **Packaging**

Please insert the cleaned and disinfected instruments in the corresponding bur blocks (if required) and pack them in single-use sterilization packagings (single or double packaging), which fulfill the following requirements (material/process):

- › EN ISO/ANSI AAMI ISO 11607 (for USA: FDA clearance)
- › suitable for steam sterilization (temperature resistance up to at least 142 °C (288 °F), sufficient steam permeability)
- › sufficient protection of the instruments as well as of the sterilization packagings to mechanical damage

##### **Sterilization**

Please use for sterilization only the listed sterilization procedures; other sterilization procedures must not be applied.

##### **Steam sterilization**

- › fractionated vacuum/dynamic air removal procedure<sup>2,3</sup> (with sufficient product drying<sup>5</sup>)
- › steam sterilizer according to EN 13060/EN 285 or ANSI AAMI ST79 (for USA: FDA clearance)
- › validated according to EN ISO 17665 (valid IQ/OQ (commissioning) and product specific performance qualification (PQ))
- › maximum sterilization temperature 138 °C (280 °F; plus tolerance according to EN ISO 17665)
- › sterilization time (exposure time at the sterilization temperature):

Area	fractionated vacuum/dynamic air removal	gravity displacement
USA	at least 4 min at 132 °C (270 °F), drying time at least 20 min <sup>4</sup>	not recommended
other countries	at least 3 min <sup>3</sup> at 132 °C (270 °F) / 134 °C (273 °F), drying time at least 20 min <sup>4</sup>	not recommended

<sup>2</sup> at least three vacuum steps

<sup>3</sup> The less effective gravity displacement procedure must not be used in case of availability of the fractionated vacuum procedure, will require significantly longer sterilization times and is to be validated dependent on product, packaging, sterilizer, program, and parameters under sole responsibility of the user.

<sup>4</sup> The effectively required drying time depends directly on parameters in sole responsibility of the user (load configuration and density, sterilizer conditions, ...) and by this is to be determined by the user. Nevertheless, drying times less than 20 min must not be applied.

<sup>5</sup> respectively 18 min (inactivation of prions, not relevant for USA)

The fundamental suitability of the instruments for an effective steam sterilization was demonstrated by an independent, governmentally accredited and recognized (§ 15 (5) MPG) test laboratory by application of the steam sterilizer HST 6x6x6 (Zirbus technology GmbH, Bad Grund) and the fractionated vacuum/dynamic air removal procedure. For this, typical conditions in clinic and doctor's practice as well as the specified procedure were considered.

The flash/immediate use sterilization procedure must not be used. Do not use dry heat sterilization, radiation sterilization, formaldehyde and ethylene oxide sterilization, as well as plasma sterilization.

##### **Storage**

Please store the instruments after sterilization in the sterilization packagings at a dry and dust-free place.

##### **Material resistance**

Please take care that the listed substances are not ingredients of the cleaning or disinfection detergent:

- › organic, mineral, and oxidizing acids (minimum admitted pH-value 5.5)
- › strong lyes (maximum admitted pH-value 11, neutral/enzymatic or alkaline cleaner recommended)
- › organic solvents (for example: acetone, ether, alcohol, benzene)
- › oxidizing agents (for example: hydrogen peroxide)
- › halogens (chlorine, iodine, bromine)
- › aromatic, halogenated hydrocarbons

<sup>6</sup>For the bur blocks alkaline cleaners must not be applied (maximum admitted pH-value 9).

Please do not clean any instruments and bur blocks by use of metal brushes or steel wool.

Please do not expose any instruments and bur blocks to temperatures higher than 142 °C (288 °F)!

Please do not apply acidic neutralizing agents or cleaning aids.

##### **Reusability**

The instruments can be reused – in case of adequate care and if they are undamaged and clean as indicated in chapter "specific aspects". The user is responsible for each further use as well as for the use of damaged and dirty instruments (no liability in case of disregard).

**Attachment A: Specific aspects**

Diamond products and ceramic grinding tools:

- › Use particular care during the cleaning of the grinding surfaces and ensure that all residues are removed

Bur blocks/instrument trays:

- › Cleaning and disinfecting only without products being loaded (products must not be cleaned and disinfected whilst they are in the bur block/ instrument tray)

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*The document in other languages as well as the latest and applicable version you find on the website www.JOTA.ch.*

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Instrument group	brush	specific/additional procedure in case of				packing	sterilization	maximum admitted cycle number (confirmed by validation, but dependent on specific application)	recommended classification according to KRINKO/RKI/BfArM guidance (only German, with respect to intended use)
		pretreat- ment	manual cleaning/ disin- fection	automated cleaning/ disinfection	maintenance				
stainless steel instruments	standard	standard	standard	standard	lubrication <u>not</u> admitted	standard	standard	10	critical B
regular steel instruments	standard	standard	standard	standard	lubrication recommended	standard	standard	10	critical B
silicone polisher	standard	standard	standard	standard	lubrication <u>not</u> admitted	standard	standard	5	critical B
endodontic instruments without stopper	endodontic brush	standard	standard	standard	lubrication <u>not</u> admitted	use of bur blocks not admitted	use of bur blocks not admitted	10	critical B
endodontic instruments with stopper	endodontic brush	mounted	mounted move the stopper at least three times during disinfection	mounted	lubrication <u>not</u> admitted	use of bur blocks not admitted	use of bur blocks not admitted	10	critical B
all other instruments	standard	standard	standard	standard	lubrication <u>not</u> admitted	standard	standard	10	critical B

**General application and safety instructions**

for the medical device from JOTA AG

Issued: January 2017

- › JOTA AG products (dental, maxillary surgery, general surgery,) must only be used by dentists, doctors and/or the respective experts who, due to their training and experience, are intensely familiar with the use of these products and who have the corresponding expertise in the respective specialist fields. The use of surgical products requires relevant expertise and experience in dental implantology, maxillary surgery and/or other surgical fields including diagnosis, preoperative planning and surgical techniques.
- › It is the sole responsibility of the doctor in charge who, depending on the respective situation (indication), decides on the actual use of the JOTA AG products for each individual case
- › All JOTA AG products have been developed for specific applications. Therefore, inappropriate use can result in the premature wear and tear of the products and put patients and users at risk.

**Application**

- › In order to avoid damaging the instruments, they must be removed from the blister pack by pulling off the back-sheet.
- › It is essential to only use turbines as well as hand and angle pieces that are technically and hygienically faultless, maintained and cleaned.
- › The instruments must be rotating when applied on material. They should not be placed on material and then brought to rotation.
- › Rotating instruments need to be clamped as far down as possible with their speed set before applying them on the object. are used with the rotary instruments.
- › Using the instruments for canting or leveraging should be avoided as it increases the risk of breakage.
- › Depending on the application, it is recommended to use protective goggles while using the instruments. Users of diamond disks should use a disk protection device.
- › Inappropriate use of the products leads to badly executed work and increased risk.
- › When working with dry materials, it is recommended to use a suction cleaning device.
- › In particular, users of hand tools should take care to use them gently and with consideration.
- › The user must at all times avoid touching the instruments and parts without protection (protective gloves should be worn).
- › Thermal bone damage caused by rotating and oscillating tools (e.g. pilot bur, conical bur, expansion bur) must at all times be avoided (user training, working at low speed and with sufficient cooling).
- › During intraoral application attention has to be made to the fact that the products are protected against aspiration or falling on the floor.

**Use of pressure**

- › Users of the instruments should at all times avoid applying excessive pressure. This can damage the working part of the instruments and cause the cutting edges to break off. At the same time, it generates excessive heat.

- › The use of excessive pressure when using grinding tools can cause the abrasive particles to break off or the instrument to become clogged and lead to heat generation.
- › During polishing, excess pressure can lead to heat generation.
- › Due to overheating, excess pressure can damage the dental pulp or, due to broken off cutting edges, it can result in undesired rough surfaces. In such cases, even instrument breakage cannot be excluded.

**Cooling**

- › In order to avoid excessive heat generation during preparation, a sterile water/sodium chloride solution supplied via a permanent external feeding device should be used to ensure sufficient cooling during use of the instruments.
- › When using FG instruments that are more than 22 mm long or whose head diameter exceeds 2 mm, additional external cooling is required.
- › Insufficient cooling will lead to irreversible damage to the bone and/or the adjacent tissue.

**Storage, disinfection, cleaning and sterilization**

- › Unless there is explicit information to the contrary, all JOTA AG products are supplied in non-sterile packaging and, depending on the application, they need to be sterilised prior to use. Prior to their first use on the patient and immediately after each use, all products need to be disinfected and sterilised. Inappropriate cleaning and sterilising of the instruments can result in the patient being infected with harmful bacteria.
- › You will find detailed instructions for the disinfecting, cleaning and sterilising of products in the Instructions for the processing of instruments produced by JOTA AG on the previous page. We would also be happy to provide you with these instructions at your request. They are also available on the internet at [www.JOTA.ch](http://www.JOTA.ch).
- › The products should be stored in appropriate, hygienically maintained containers. The same applies to sterilised instruments. The stored products must be protected from dust, humidity and recontamination. Instructions as to maximum storage duration must be adhered to.

**Speed recommendations for rotary instruments**

- › Following the instrument-specific speed recommendations produces the best results.
- › Exceeding the maximum admissible speed (rpm) when using long and pointed instruments tends to produce vibrations that can lead to the destruction of the instrument.
- › When using working parts with diameters exceeding the thickness of the shaft, excessive speed can release great centrifugal forces that may cause the shaft to bend and/or the instrument to break. Therefore, the maximum admissible rpm must never be exceeded.
- › Please consult the manufacturer's information (see catalogue or [www.JOTA.ch](http://www.JOTA.ch)) for the recommended and the maximum admissible speed ranges. Non-compliance with the maximum admissible speed puts safety at risk.
- › Generally, the following rules apply:

- The larger the working part of an instrument the lower the speed
- Surgical instruments: suitable for geared down micro-motor hand and angle pieces 10:1 with stable ball bearings. Speed 600 to 800 rpm with physical and, possibly, sterile external cooling or internal cooling when using the respective hand piece.

**Discarding worn instruments and parts**

- › JOTA AG products can principally be reused several times – unless specifically indicated and labelled otherwise. Rotating instruments are subject to wear. The option of and accountability for multiple use of product and the frequency of application is solely the decision and own responsibility of the treating clinician based on the application in each case and the possible wear of the products. If in doubt, the products should always be sorted out early and replaced.
- › Broken off cutting edges of instruments cause vibrations and great forces of pressure, which, in turn, leads to broken preparation corners and rough surfaces.
- › Bare patches on diamond instruments indicate a lack of abrasive particles and can be a sign of blunt instruments. This leads to excessive temperatures during instrument use.
- › Instruments that are bent and/or do not run true should be discarded forthwith.
- › With the reuse of disposable products the risk of infection cannot be excluded and a risk-free functional safety cannot be guaranteed.

**Additional instructions regarding the use of trepans**

- › When using trepans, you have to proceed with particular care. For example, it is advisable not to exceed the recommended rpm speed ranges.
- › In order to prepare for the actual use of a trepan, it should be set to produce counter-clockwise rotations creating a groove in the bone. Afterwards the trepan can be inserted into this groove and, using clockwise rotations, it can be moved further down.
- › Carrying out a prior X-Ray is essential to establish the maximum possible drilling depth and to maintain the necessary distance, for example, to the mandibular nerve. As an additional safety measure to spare the nerve, the axial direction of the trepan countersink attachment, based on the sagittal level of the ascending branch, must be milled laterally at an angle of approx. 15-20°.

$\varnothing < 045$

$\varnothing \geq 045$

**Further comments:**

- › Due to statutory regulations, returned goods can, on principle, only be accepted if the complete batch number is provided. This number can be found on the product packaging.

jota

**JOTA AG** Rotary Instruments

Hirschenprungstrasse 2, 9464 Rüthi, Switzerland  
Phone +41 (0)71 767 79 99, Fax +41 (0)71 767 79 97  
info@JOTA.ch, [www.JOTA.ch](http://www.JOTA.ch) + SWITZERLAND